

ThermoFisher ND-1000 NanoDrop

1. Overview:

- a. The Thermo Scientific NanoDrop™ 1000 Spectrophotometer measures 1 ul samples with high accuracy and reproducibility. The full spectrum (220nm-750nm) spectrophotometer utilizes a patented sample retention technology that employs surface tension alone to hold the sample in place.
- b. The small sample requirement and ease of use make the NanoDrop 1000 Spectrophotometer ideally suited for measuring:
 - Nucleic acid concentration and purity of nucleic acid samples up to 3700 ng/ul (dsDNA) without dilution
 - Fluorescent dye labeling density of nucleic acid microarray samples
 - Purified protein analysis (A280) up to 100 mg/ml (BSA)
 - Expanded spectrum measurement and quantitation of fluorescent dye labeled proteins, conjugates, and metalloproteins
 - Bradford Assay analysis of protein
 - BCA Assay analysis of protein
 - Lowry Assay analysis of protein
 - Pierce Protein 660 nm Protein Assay
 - Cell density measurements
 - General UV-Vis spectrophotometry

2. Training and Use:

- a. Before using the NanoDrop, new users must be trained. Contact admin.tbep@utoronto.ca to arrange a training session.
- b. Book your timeslot in Google calendar under the account bookhere.equipment@gmail.com, and indicate your name and the lab to which you belong.

3. Responsibility:

- a. Labs are responsible for providing their own consumables (gloves, Kimwipes, etc.) for use with the NanoDrop.
- b. Users must properly clean the NanoDrop after use.

4. Precautions:

- a. **ALWAYS** wear gloves when touching the instrument.
- b. Keep stainless steel arm closed whenever possible; avoid dust contamination.
- c. Do **NOT** place a kimwipe between the pedestals when not using the Nanodrop.

5. Procedure:

a. Preparation:

- Quickly wipe the upper and lower pedestal with a Kimwipe soaked in RNaseZAP.
- Promptly rinse off the RNaseZap from the pedestals with a Kimwipe that has been soaked in water from the "Cleaning" bottle (Ultrapure or RNase/DNase-free dH₂O).

b. Initialization:

- Choose "Nucleic Acid". Initialize the instrument by following the on-screen prompt to load a water sample (1.5-2 uL) from the "Blanking" water bottle (UltraPure or RNase/DNase-free dH₂O). Click "OK".
- Once the initialization is complete, vigorously wipe the upper and lower pedestal dry with a clean Kimwipe. Between all subsequent measurements, repeat this cleaning.

c. Quantification:

- For RNA quantification, choose "RNA-40".
- Load a 1.5-2 uL buffer blank (typically water, you can use the "Blanking" water bottle). Click "Blank".
- Load a 1.5-2 uL buffer blank, click "Measure".
- Confirm that readings near 0 ng/uL are measured.

d. Measuring:

- Load 1.5-2 uL of your first sample. Enter a sample ID. Click "Measure".
- Repeat for all your samples, wiping the pedestal clean every time.

e. Data Export:

- Export your data by clicking "Show Report">"Reports" (top toolbar)>"Save Report...">"Export Report & Standards Tables". Save the report in the following directory: C:/nanodrop data/Reports/your username folder.

f. Cleaning and Maintenance:

- After you have measured your final sample, wipe the upper and lower pedestal dry with a clean Kimwipe. Apply 5 ul of "Cleaning" water onto each bottom pedestal.
- Lower the upper pedestal arm to form a liquid column; let it sit for approximately 2-3 minutes.
- Wipe both pedestals dry one more time with a new Kimwipe.
- NOTE: Typically dH₂O is sufficient for removal of samples that have dried on the optical pedestals. There are a few cases (i.e. dried proteins) that may require a more rigorous cleaning protocol. For these cases, we recommend that 0.5M HCl be substituted for the 5 ul of dH₂O. Please follow an HCl application with 5 ul of dH₂O to ensure any residual HCl is removed.

- **REMINDER:** Do **NOT** place a kimwipe between the pedestals when not using the Nanodrop. Leave the pedestal **CLOSED** after use.

6. References:

- a. Refer to the ND-1000 User Manual for additional information.